



DOBB[®] **Sampling Device**

DOBB[®]'s simplify and improve the aim-orientated sampling of different materials (water and biota). Due to the use of high quality materials *DOBB*[®]'s are particularly suitable to collect chemical and microbiological samples.

Technical data

Our *DOBB*[®]'s, being made of high quality borosilicate glass 3.3, are available in two sizes:

***DOBB*[®] GL 45**

- Length: 200 mm
- Opening diameter (inside): 35 mm
- Sample volume (inside): 300 ml
- Weight (empty): 300 g

***DOBB*[®] GL 80**

- Length: 200 mm
- Opening diameter (inside): 60 mm
- Sample volume (inside): 800 ml
- Weight (empty): 500 g



Accessories

- Various caps and pouring rings
- Transport solution for surface and underwater
- Special sizes available on request

