

LiSSa

Sediment Sampler

LiSSa enables diving scientists to recover nearly undisturbed sediment cores by an easy and proven method.

Technical data

Available lengths of the sample cylinder:


- 100 cm
- 200 cm
- Custom lengths available on demand

Available nominal diameter (DN):

- 32 mm
- 50 mm

Material:

- PVC (sampling cylinder)
- Core catcher steel or PVC



Accessories

- Caps DN 32 and 50
- Core catcher (DN 32, DN 50) for various sample textures and grain sizes
- Lift bag to recover lengthy liners and those that are stuck in adhesive substrate

